



Catalog Number: 10890-R211

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human BAMBI / NMA Protein (Catalog#10890-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human BAMBI / NMA (rh BAMBI / NMA; Catalog#10890-H08H; NP_036474.1; Met 1-Ala 152).
Ig Type:	Rabbit IgG
Clone ID:	211
Specificity:	Human BAMBI / NMA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human BAMBI.

Please Note: Optimal concentrations/dilutions should be determined by the end user.