

CCL5 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10900-MM42

GENERAL INFORMATION

Immunogen:	Recombinant Human CCL5 Protein (Catalog#10900-HNAE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CCL5 (rh CCL5; Catalog#10900-HNAE; NP_002976.2; Ser24-Ser91). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	42
Specificity:	Human CCL5
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	D17S136E, eoCP, RANTES, SCYA5, SIS-delta, SISd, TCP228

APPLICATIONS

Applications:	ELISA(Cap)
----------------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human CCL5 in ELISA pair set(Catalog: # SEK10900). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 10900-MM05).
------------------------------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.