

FGF-21 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10911-MM25

GENERAL INFORMATION

Immunogen:	Recombinant Human FGF-21 Protein (Catalog#10911-H07E)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human FGF-21 (rh FGF-21; Catalog#10911-H07E; NP_061986.1; His29-Ser209). The IgG fraction of t
Ig Type:	Mouse IgG2a
Specificity:	Human FGF-21
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000
-------	----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.