



Catalog Number: 10924-R009


EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human FSTL1 protein (Catalog#10924-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human FSTL1 (rh FSTL1; Catalog#10924-H08H; NP_009016.1; Met1-Ile308).
Ig Type:	Rabbit IgG
Clone ID:	009
Specificity:	Human FSTL1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	FRP,FSL1,MIR198
APPLICATIONS	
Applications:	FCM,ICC/IF
RECOMMENDED CONCENTRATION	
ICC/IF	ICC/IF: 1:20-1:100
Flow Cytometry	FCM: 1:25-1:100


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Catalog Number: 10924-R009

EliteRmab® is a registered trademark of Sino Biological Inc.



Confocal immunofluorescence staining of Human FSTL1 in HeLa cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human FSTL1 monoclonal antibody (1:60) at 37°C 1 hour. Then cells were stained with the Alexa Fluor® 488-conjugated goat Anti-rabbit IgG secondary antibody (green), counterstained with DAPI (blue). Positive staining was localized to cytoplasm.



Flow cytometric analysis of Human FSTL1 expression on HeLa cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with purified anti-Human FSTL1, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.