

GDF15/MIC-1 Antibody (HRP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10936-MM12-H

GENERAL INFORMATION

Immunogen:	Recombinant Human GDF15/MIC-1 Protein (Catalog#10936-H07Y)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human GDF15/MIC-1 (rh GDF15/MIC-1; Catalog#10936-H07Y; NP_004855.2; Ala197-Ile308). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Mouse IgG1
Clone ID:	12
Specificity:	Human GDF15/MIC-1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	GDF15

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.1-1 µg/ml
-------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>