

# CD2AP Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 109645-T32

## GENERAL INFORMATION

|                     |                                                                                                                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen:</b>   | A synthetic peptide corresponding to the C-terminus of the Human CD2AP                                                                                                                                                                             |
| <b>Preparation</b>  | Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human CD2AP, and purified by antigen affinity chromatography.                                                                                        |
| <b>Ig Type:</b>     | Rabbit IgG                                                                                                                                                                                                                                         |
| <b>Specificity:</b> | Human CD2AP                                                                                                                                                                                                                                        |
| <b>Formulation:</b> | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                    |
| <b>Storage:</b>     | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |

## APPLICATIONS

|                      |    |
|----------------------|----|
| <b>Applications:</b> | WB |
|----------------------|----|

## RECOMMENDED CONCENTRATION

|                     |                  |
|---------------------|------------------|
| <b>Western Blot</b> | WB: 1:500-1:2000 |
|---------------------|------------------|


**Please Note:** Optimal concentrations/dilutions should be determined by the end user.

# CD2AP Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 109645-T32



Anti-CD2AP rabbit polyclonal antibody at 1:500 dilution

Lane A: HeLa Whole Cell Lysate

Lane B: Jurkat Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 71 kDa

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

**Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)**

**Website: <http://www.sinobiological.com>**