

CD8a / Lyt2 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10980-MM22

GENERAL INFORMATION

Immunogen:	Recombinant Human CD8a / Lyt2 protein (Catalog#10980-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD8a / Lyt2 (rh CD8a / Lyt2; Catalog#10980-H08H; NP_001759.3; Met1-Asp182). The IgG fraction of the cell culture supernatant was purified by Protein A affinity.
Ig Type:	Mouse IgG1
Clone ID:	22
Specificity:	Human CD8a / Lyt2
	No cross-reactivity in ELISA with Human cell lysate (293 cell line)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human CD8a / Lyt2.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.