

CD20/MS4A1 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11007-MM06

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD20/MS4A1 Protein (Catalog#11007-H07H2)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD20/MS4A1 (rh CD20/MS4A1; Catalog#11007-H07H2; NP_068769.2; Glu213-Pro297). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	06
Specificity:	Human, Cynomolgus
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:30-1:100


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD20/MS4A1 Antibody, Mouse MAb




Sino Biological
Biological Solution Specialist


Catalog Number: 11007-MM06




Immunohistochemical staining CD20 in human tonsil with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining CD20 in human lymphoma with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining CD20 in cynomolgus spleen with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining CD20 in cynomolgus lymph node with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.