

CD20 / MS4A1 Antibody (PerCP), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11007-R001-C

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Human CD20 Transfected CHO Cells
Reagents:	PerCP-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, Human CD20 Transfected CHO Cells and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.
Ig Type:	Rabbit IgG
Clone ID:	001
Specificity:	Human CD20 / MS4A1
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	MS4A1

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.


CD20 / MS4A1 Antibody (PerCP), Rabbit MAb

Catalog Number: 11007-R001-C



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Profile of anti-CD20 reactivity on peripheral blood lymphocytes analyzed by flow cytometry.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.