


Catalog Number: 110122-R060


EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	A synthetic peptide corresponding to the Human GCG/Glucagon
Preparation	This antibody was obtained from a rabbit immunized with purified, synthetic peptide corresponding to the Human GCG/Glucagon.
Ig Type:	Rabbit IgG
Clone ID:	060
Specificity:	Human GCG/Glucagon
	Reacts with: Human, Mouse, Rat Pig, Dog, Bovine (Species predicted to react based on 100% sequence homology)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, ELISA(Det), IHC-P
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000
IHC-P	IHC-P: 1:1000-1:5000
ELISA(Det)	ELISA(Det): 1:1000-1:10000 In a sandwich ELISA, one of antibody (Cat# 110122-R035, 110122-R060) can be used as detection antibody when paired with Cat# 110122-R030.


Please Note: Optimal concentrations/dilutions should be determined by the end user.



Immunohistochemical staining of human Glucagon in human pancreas with rabbit monoclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of Mouse Glucagon in Mouse pancreas with rabbit monoclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of rat cagon in rat pancreas with rabbit monoclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>