



Catalog Number: 110182-R014

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the Mouse Ins1/Insulin-1
Preparation	This antibody was obtained from a rabbit immunized with purified, synthetic peptide corresponding to the Mouse Ins1/Insulin-1
Ig Type:	Rabbit IgG
Clone ID:	014
Specificity:	Mouse Ins1/Insulin-1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA, ELISA(Det)
----------------------	-------------------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
ELISA(Det)	ELISA(Det): 1:1000-1:10000 In a sandwich ELISA, one of antibody (Cat# 110182-R012, 110182-R014) can be used as detection antibody when paired with Cat# 110182-R029.

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>