

# CD8B / P37 / LEU2 Antibody, Rabbit MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 11031-R025

EliteRmab® is a registered trademark of Sino Biological Inc.

| GENERAL INFORMATION       |                                                                                                                                                                                                                                                                                                                                                                                                     |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:                | Recombinant Human CD8B / P37 / LEU2 Protein (Catalog#11031-HCCH)                                                                                                                                                                                                                                                                                                                                    |
| Preparation               | This antibody was obtained from a rabbit immunized with purified, recombinant Human CD8B / P37 / LEU2 (rh CD8B / P37 / LEU2; Catalog#11031-HCCH; P10966-1; Met1-Pro170).                                                                                                                                                                                                                            |
| Ig Type:                  | Rabbit IgG                                                                                                                                                                                                                                                                                                                                                                                          |
| Clone ID:                 | 025                                                                                                                                                                                                                                                                                                                                                                                                 |
| Specificity:              | Human CD8B / P37 / LEU2                                                                                                                                                                                                                                                                                                                                                                             |
| Formulation:              | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                                                                                                                                                                     |
| Storage:                  | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.<br>Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles. |
| APPLICATIONS              |                                                                                                                                                                                                                                                                                                                                                                                                     |
| Applications:             | ELISA                                                                                                                                                                                                                                                                                                                                                                                               |
| RECOMMENDED CONCENTRATION |                                                                                                                                                                                                                                                                                                                                                                                                     |
| ELISA                     | ELISA: 1:5000-1:10000<br>This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CD8B / P37 / LEU2.                                                                                                                                                                                                                                                     |

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**