

CD8B / P37 / LEU2 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11031-R039

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human CD8B / P37 / LEU2 Protein (Catalog#11031-HCCH)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD8B / P37 / LEU2 (rh CD8B / P37 / LEU2; Catalog#11031-HCCH; P10966-1; Met1-Pro170).
Ig Type:	Rabbit IgG
Clone ID:	039
Specificity:	Human CD8B / P37 / LEU2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
--------------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.


CD8B / P37 / LEU2 Antibody, Rabbit MAb

Catalog Number: 11031-R039




Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunochemical staining of human CD8B in human spleen with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).



Immunochemical staining of human CD8B in human lymph node with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).