

Insulin / INS Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11038-T08

GENERAL INFORMATION

Immunogen:	Recombinant Human Insulin / INS Protein (Catalog#11038-HNAY)
Preparation	Produced in rabbits immunized with purified, recombinant Human Insulin / INS (rh Insulin / INS; Catalog#11038-HNAY; P01308; Phe25-Lys53 & Gly90-Asn110). Insulin / INS specific IgG was purified by Human Insulin / INS affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Insulin / INS
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
--------------	---------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Insulin / INS Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11038-T08



Immunohistochemical staining of human Insulin in human pancreas with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).