

**Influenza A H1N1 (A/California/04/2009)
Hemagglutinin / HA Antibody (HRP),
Rabbit PAb**

Catalog Number: 11055-RP01-H



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA Protein
Preparation	Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA. Total IgG was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Rabbit IgG
Specificity:	Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Hemagglutinin,HA

APPLICATIONS

Applications:	ELISA
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.1-1 µg/mL
--------------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.