

Influenza A H3N2 Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 11056-T62



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION	
Immunogen:	Recombinant H3N2 (A/Brisbane/10/2007) HA protein (Catalog#11056-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant H3N2 (A/Brisbane/10/2007) HA / Hemagglutinin (Catalog#11056-V08H; ABW23353.1; Met 1-Ile 531). Total IgG was purified by Protein A affinity chromatography.
Specificity:	Influenza A H3N2 (A/Brisbane/10/2007) HA protein
	Has cross-reactivity in WB and ELISA with H3N2 (A/Darwin/1415/2025) HA trimer Protein (Cat# 41128-V08H) Has cross-reactivity in ELISA with H1N1 (A/Wisconsin/67/2022) HA Trimer Protein (Cat# 40940-V08H2) H1N1 (A/Wisconsin/67/2022) HA Protein (Cat# 40940-V08B) H1N1 (A/Wisconsin/67/2022) HA Protein (Cat# 40940-V08H) H3N2 (A/Darwin/9/2021) HA Protein (Cat# 40859-V08B) H3N2 (A/Darwin/9/2021) HA Protein (Cat# 40859-V08H) H3N2 (A/Massachusetts/18/2022; A/Thailand/8/2022) HA Protein (Cat# 40992-V08B) H3N2 (A/Massachusetts/18/2022; A/Thailand/8/2022) HA Protein (Cat# 40992-V08H) Influenza B (B/Austria/1359417/2021) HA Protein (Cat# 40862-V08H) Influenza B (B/Austria/1359417/2021) HA Protein (Cat# 40862-V08B)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2 °C -8 °C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20 °C to -80 °C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	HA
APPLICATIONS	
Applications:	WB, ELISA
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:1000-1:5000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Influenza A H3N2 Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 11056-T62



Sino Biological
Biological Solution Specialist




Anti-HA rabbit polyclonal antibody at 1:5000 dilution.

Sample: Influenza A H3N2
(A/Brisbane/10/2007) Hemagglutinin /
HA Recombinant Protein
Lane A: 30ng
Lane B: 5ng

Secondary
Goat Anti- Rabbit IgG H&L (Dylight
800) at 1/10000 dilution.

Developed using the ECL technique.
Performed under reducing conditions.




Anti-HA rabbit polyclonal antibody at 1:1000 dilution.

Sample: Recombinant Protein 10 ng
Lane 1: A/Croatia/10136RV/2023 (H3N2) HA
(Cat#41058-V08B)
Lane 2: A/Croatia/10136RV/2023 (H3N2) HA
(Cat#41058-V08H)
Lane 3: A/District of Columbia/27/2023 (H3N2)
HA (Cat#41056-V08H)
Lane 4: A/District of Columbia/27/2023 (H3N2)
HA (Cat#41056-V08B)

Secondary
Goat Anti-rabbit IgG (H+L)/HRP at 1/10000
dilution

Developed using the ECL technique.
Performed under reducing conditions



Anti-Influenza A H3N2 (A/Brisbane/10/2007)
Hemagglutinin/HA rabbit polyclonal antibody
at 1:1000 dilution

Lane A: Influenza A H3N2 (A/New
York/GKISBBBG87773/2025) Hemagglutinin
/ HA trimer Protein(Cat: 41128-V08H) 10ng

Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.
Performed under reducing conditions.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>