

Phospho-Met (Tyr1234, 1235) Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 110587-R0060

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to residues around (Tyr1234, 1235) of Human Phospho-Met
Preparation	This antibody was obtained from a rabbit immunized with a synthetic peptide corresponding to residues around (Tyr1234, 1235) of Human Phospho-Met.
Ig Type:	Rabbit IgG
Clone ID:	0060
Specificity:	Human Species predicted to react based on 100% sequence homology: Mouse, Rat, Cynomolgus, Rabbit
Formulation:	10 mM sodium HEPES, 150 mM NaCl, 100 µg/mL BSA, 50% glycerol, pH 7.5
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB
----------------------	----

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:1000-1:5000
---------------------	-------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.


Phospho-Met (Tyr1234, 1235) Antibody, Rabbit MAb

Catalog Number: 110587-R0060




Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Western blot analysis of extracts from serum-starved A549, untreated (line A); treated with HGF (Cat# 10463-HNAS) (50 ng/mL, 5min; +) (line B); treated with HGF and λ -phosphatase (Cat# L01-20U) (line C) using Phospho-Met (Tyr1234, 1235) rabbit monoclonal Antibody (Cat# 110587-R0060) at 1:1000 dilution.

(Validation Experiment)



Western blot analysis of extracts from A549, untreated (line A) or treated with HGF (Cat# 10463-HNAS) (50 ng/mL, 5min; +) (line B), using Phospho-Met (Tyr1234, 1235) rabbit monoclonal Antibody (Cat# 110587-R0060) at 1:1000 dilution (upper) or Anti-Actin Antibody, Rabbit Monoclonal (Cat# 100166-MCR10) at 1:50000 dilution (lower).