

IL32 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11064-R403

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human IL32 Protein (Catalog#11064-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL32 (rh IL32; Catalog#11064-H08H; NP_001012651.1; Met1-Lys131).
Ig Type:	Rabbit IgG
Clone ID:	403
Specificity:	Human IL32
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA(Det)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human IL32 in ELISA pair set (Catalog: # SEK11064). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 11064-R034).

Please Note: Optimal concentrations/dilutions should be determined by the end user.