

Coagulation Factor X/F10 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11076-MM10

GENERAL INFORMATION

Immunogen:	Recombinant Human Coagulation Factor X/F10 Protein (Catalog#11076-H08B)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Coagulation Factor X/F10 (rh Coagulation Factor X/F10; Catalog#11076-H08B; NP_000495.1; Met1-Lys488). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Clone ID:	10
Ig type:	Mouse IgG1
Specificity:	Human Coagulation Factor X/F10
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA(Cap)
----------------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human Coagulation Factor X/F10 in ELISA pair set (Catalog: # SEK11076). it can be used as capture antibody when paired with (Catalog: # 11076-R058).
--	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.