

CEACAM5 Antibody (HRP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11077-MM07-H

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD66e protein (Catalog#11077-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CEACAM5 / CD66e (rh CEACAM5 / CD66e; Catalog#11077-H08H; NP_004354.2; Met 1-Ala 685). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Mouse IgG1
Clone ID:	1B8E7
Specificity:	Human CEACAM5
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.