

IL18R1 Antibody (HRP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11102-MM05-H

GENERAL INFORMATION

Immunogen:	Recombinant Human IL18R1 Protein (Catalog#11102-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IL18R1 (rh IL18R1; Catalog#11102-H08H; NP_003846.1; Met1-Arg329). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Mouse IgG1
Clone ID:	05
Specificity:	Human IL18R1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	IL18R1

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.1-1 µg/ml
-------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.