

IL18R1 / CD218a Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11102-MM08

GENERAL INFORMATION	
Immunogen:	Recombinant Human IL18R1 protein (Catalog#11102-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IL18R1 / CD218a (rh IL18R1; Catalog#11102-H08H; NP_003846.1; Met 1-Arg 329). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	9A8D6F3
Specificity:	Human IL18R1 / CD218a
	No cross-reactivity in ELISA with Human IL1R1 / CD121a Human IL1R2 / CD121b Human IL1R3 / IL1RAP Human IL1RL1 / IL1R4 / ST2 Human IL18RAP / IL1R7 Human IL1RAPL1 / IL1R8 Human IL1R9 / IL1RAPL2 Mouse IL18R1 / CD218a Mouse TIR8 / SIGIRR Human cell lysate (293 cell line)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human IL18R1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.