

Phospho-ATR (Ser428) Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 111263-R0025

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to residues around Ser428 of Human ATR
Preparation:	This antibody was obtained from a rabbit immunized with a synthetic peptide corresponding to residues around Ser428 of Human ATR.
Ig Type:	Rabbit IgG
Clone ID:	0025
Specificity:	Human Phospho-ATR (Ser428)
Formulation:	10 mM sodium HEPES, 150 mM NaCl, 100 µg/mL BSA, 50% glycerol, pH 7.5
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB
----------------------	----

RECOMMENDED CONCENTRATION



Western Blot	WB: 1:500-1:2000
---------------------	------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>



(Validation Experiment)