

Phospho-cdc25C (Ser216) Antibody, Rabbit MAb




Sino Biological
Biological Solution Specialist

Catalog Number: 111281-R0043


EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	A synthetic phosphopeptide corresponding to residues around Ser216 of the Human Phospho-cdc25C
Preparation	This antibody was obtained from a rabbit immunized with a synthetic phosphopeptide corresponding to residues around Ser216 of the Human Phospho-cdc25C.
Ig Type:	Rabbit IgG
Clone ID:	0043
Specificity:	Human Phospho-cdc25C (Ser216)
	Reacts with: Human
Formulation:	10 mM sodium HEPES, 150 mM NaCl, 100 µg/mL BSA, 50% glycerol, pH 7.5
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1000-1:5000

Please Note: Optimal concentrations/dilutions should be determined by the end user.



Western blot analysis of extracts from HT-29(-) and serum-starved HT-29 (+) using Phospho-cdc25C (Ser216) rabbit monoclonal Antibody (Cat# 111281-R0043) at 1:1000 dilution (upper) or Beta-Tubulin Loading Control Antibody, Mouse Mab (Cat: 100109-MM05T) at 1:10000 dilution (lower).



Western blot analysis of extracts from HT-29 (line A); serum-starved HT-29 (line B); serum-starved HT-29 treated with λ-phosphatase (line C) using Phospho-cdc25C (Ser216) rabbit monoclonal Antibody (Cat# 111281-R0043) at 1:1000 dilution.

(Validation Experiment)

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>