

B7-H3 / CD276 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11188-MM06-P

GENERAL INFORMATION

Immunogen:	Recombinant Human B7-H3 / CD276 protein (Catalog#11188-H08H)
Reagents:	PE-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human B7-H3 / CD276 (rh B7-H3 / CD276; Catalog#11188-H08H; Q5ZPR3-1; Met 1-Thr 461) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG2b
Clone ID:	06
Specificity:	Human B7-H3 / CD276
Concentration:	2 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	CD276

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

B7-H3 / CD276 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11188-MM06-P



Flow cytometric analysis of anti-B7H3 (CD276) reactivity on PC-3 cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.