

B7-H3/CD276 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11188-T24

GENERAL INFORMATION	
Immunogen:	Recombinant Human B7-H3/CD276 Protein (Catalog#11188-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human B7-H3/CD276 (rh B7-H3/CD276; Catalog#11188-H08H; Q5ZPR3-1; Met1-Thr461). B7-H3/CD276 specific IgG was purified by Human B7-H3/CD276 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human B7-H3/CD276
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:2000-1:10000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human B7-H3/CD276.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

B7-H3/CD276 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 11188-T24



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human B7H3 in human ovarian cancer with rabbit polyclonal antibody at 1:5000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human B7H3 in human placenta with rabbit polyclonal antibody at 1:5000 dilution, formalin-fixed paraffin embedded sections.