

Phospho-Stat2 (Tyr690) Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 111963-R0001

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to residues around Tyr690 of Human Stat2
Preparation:	This antibody was obtained from a rabbit immunized with a synthetic peptide corresponding to residues around Tyr690 of Human Stat2.
Ig Type:	Rabbit IgG
Clone ID:	R0001
Specificity:	Human Species predicted to react based on 100% sequence homology: Cynomolgus Phospho-Stat2 (Tyr690)
Formulation:	10 mM sodium HEPES, 150 mM NaCl, 100 µg/mL BSA, 50% glycerol, pH 7.5
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB
----------------------	----

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:10000-1:100000
---------------------	----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>


Phospho-Stat2 (Tyr690) Antibody, Rabbit MAb




Sino Biological
Biological Solution Specialist

Catalog Number: 111963-R0001

EliteRmab® is a registered trademark of Sino Biological Inc.



Western blot analysis of extracts from serum-starved HeLa, untreated (line A) or treated with IFN α (13833-HNAY) (10ng/ml, 30min; +) (line B), using Phospho-Stat2 (Tyr690) rabbit monoclonal Antibody (Cat: 111963-R0001) at 1:20000 dilution.



Western blot analysis of extracts from serum-starved HeLa, untreated (line A) or treated with IFN α (13833-HNAY) (10ng/ml, 30min; +) (line B); treated with IFN α and λ -phosphatase (Cat: L01-20U) (line C) using Phospho-Stat2 (Tyr690) rabbit monoclonal Antibody (Cat: 111963-R0001) at 1:20000 dilution.
(Validation Experiment)

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>