

Apolipoprotein H Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11221-MM08

GENERAL INFORMATION

Immunogen:	Recombinant Human Apolipoprotein H Protein (Catalog#11221-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Apolipoprotein H (rh Apolipoprotein H; Catalog#11221-H08H; NP_000033.2; Met1-Cys345). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	08
Specificity:	Human Apolipoprotein H
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA(Det)
---------------	------------

RECOMMENDED CONCENTRATION

ELISA(Det)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human Apolipoprotein H in ELISA pair set(Catalog: # SEK11221). In a sandwich ELISA, it can be used as detection antibody.
------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.