

GPA33 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11277-RP02

GENERAL INFORMATION	
Immunogen:	Recombinant Human GPA33 protein (Catalog#11277-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human GPA33 (rh GPA33; Catalog#11277-H08H; Q99795-1; Met 1-Val 235). GPA33 specific IgG was purified by Human GPA33 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human GPA33
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human GPA33.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

GPA33 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 11277-RP02




Immunohistochemical staining of human GPA33 in human rectal cancer with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections). Positive staining was localized to epithelium.



Immunohistochemical staining of human GPA33 in human colon with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections). Positive staining was localized to epithelium.



Immunohistochemical staining of human GPA33 in human colon carcinoma with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections). Positive staining was localized to epithelium.



Immunohistochemical staining of human GPA33 in human appendix with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections). Positive staining was localized to epithelium.