

Adenylosuccinate Lyase Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11287-MM05T

GENERAL INFORMATION

Immunogen:	Recombinant Human Adenylosuccinate Lyase Protein (Catalog#11287-H07E)
Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells.
Ig Type:	Mouse IgG1
Clone ID:	05
Specificity:	Human Adenylosuccinate Lyase
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ICC/IF
---------------	--------

RECOMMENDED CONCENTRATION

ICC/IF	ICC/IF: 1:20-1:100
--------	--------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Adenylosuccinate Lyase Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11287-MM05T



Immunofluorescence staining of His-ADSL in HeLa cells. Cells were fixed with 4% PFA, permeabilized with 0.1% Triton X-100 in PBS, blocked with 10% serum, and incubated with mouse anti-Human His-ADSL monoclonal antibody (dilution ratio 1:60) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-mouse IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to Cytoplasm.