

Epha4 / HEK8 / TYRO1 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11314-RP01

GENERAL INFORMATION	
Immunogen:	11314-H08H
Preparation	Produced in rabbits immunized with purified, recombinant Human EPHA4 (rh EPHA4; Catalog#11314-H08H; NP_004429.1; Met 1-Thr 547). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human EPHA4 / HEK8 / TYRO1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	HEK8
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000

Please Note: Optimal concentrations/dilutions should be determined by the end user.