

BLMH / Bleomycin Hydrolase Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11333-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Human BLMH / Bleomycin Hydrolase protein (Catalog#11333-H07E)
Preparation	Produced in rabbits immunized with purified, recombinant Human BLMH / Bleomycin Hydrolase (rh BLMH; Catalog#11333-H07E; NP_000377.1; Ser2-Glu455). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human BLMH / Bleomycin Hydrolase
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human BLMH / Bleomycin Hydrolase.

Please Note: Optimal concentrations/dilutions should be determined by the end user.