

Cadherin 17 / CDH17 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11360-T26

GENERAL INFORMATION

Immunogen:	Recombinant Human Cadherin 17 / CDH17 Protein (Catalog#11360-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human Cadherin 17 / CDH17 (rh Cadherin 17 / CDH17; Catalog#11360-H08H; Q12864-1; Met1-Met787). Cadherin 17 / CDH17 specific IgG was purified by Human Cadherin 17 / CDH17 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Cadherin 17 / CDH17
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	CDH17

APPLICATIONS

Applications:	ELISA, IHC-P, ICC/IF
----------------------	----------------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:2000-1:20000
ICC/IF	ICC/IF: 1:50-1:200
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human Cadherin 17 / CDH17.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Cadherin 17 / CDH17 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 11360-T26



Immunohistochemical staining of human CDH17 in human rectum carcinoma with rabbit polyclonal antibody at 1:10000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CDH17 in human colon carcinoma with rabbit polyclonal antibody at 1:10000 dilution, formalin-fixed paraffin embedded sections.



Immunofluorescence staining of CDH17 in CACO2 cells. Cells were fixed with 4% PFA, blocked with 10% serum, and incubated with rabbit anti-Human CDH17 polyclonal antibody (dilution ratio 1:100) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to Cell membrane.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>