

**COCH Antibody, Rabbit PAb,
Antigen Affinity Purified**



Sino Biological
Biological Solution Specialist

Catalog Number: 11368-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human COCH Protein (Catalog#11368-H07H)
Preparation	Produced in rabbits immunized with purified, recombinant Human COCH (rh COCH; Catalog#11368-H07H; NP_001128530.1; Glu25-Gln550). COCH specific IgG was purified by Human COCH affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human COCH
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
--------------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.