



Catalog Number: 11389-R017

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human FLRT1 protein (Catalog#11389-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human FLRT1 (rh FLRT1; Catalog#11389-H08H; Q9NZU1-1; Met1-Pro524).
Ig Type:	Rabbit IgG
Clone ID:	017
Specificity:	Human FLRT1
	No cross-reactivity in ELISA with Human FLRT3 Human FLRT2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human FLRT1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.