

IL21 Receptor Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11483-R112

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human IL21 Receptor Protein (Catalog#11483-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL21 Receptor (rh IL21 Receptor; Catalog#11483-H08H; Q9HBE5-1; Met1-Pro236).
Ig Type:	Rabbit IgG
Clone ID:	112
Specificity:	Human IL21 Receptor
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>