

CD200RLa Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11576-RP02

GENERAL INFORMATION	
Immunogen:	Recombinant Cynomolgus CD200RLa protein (Catalog#11576-C08H)
Preparation	Produced in rabbits immunized with purified, recombinant Cynomolgus CD200RLa (rC CD200RLa; Catalog#11576-C08H; Met 1-Leu 240). CD200RLa specific IgG was purified by cynomolgus CD200RLa affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Cynomolgus CD200RLa
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Cynomolgus CD200rla.

Please Note: Optimal concentrations/dilutions should be determined by the end user.