

BirA / Bifunctional protein BirA Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11582-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant E. coli BirA / Bifunctional protein BirA protein (Catalog#11582-E10E)
Preparation	Produced in rabbits immunized with purified, recombinant E. coli BirA / Bifunctional protein BirA (Catalog#11582-E10E; P06709; Met1-Lys321). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	E. coli BirA / Bifunctional protein BirA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
	(Antibody's applications have not been validated with corresponding E. coli. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect E. coli BirA / Bifunctional protein BirA.

Please Note: Optimal concentrations/dilutions should be determined by the end user.