

GBA3 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11593-MM05-P

GENERAL INFORMATION

Immunogen:	Recombinant Human GBA3 protein (Catalog#11593-H08B)
Reagents:	PE-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human GBA3 (rh GBA3; Catalog#11593-H08B; NP_066024.1; Met1-Leu469) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG2a
Clone ID:	05
Specificity:	Human GBA3
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	GBA3

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

GBA3 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11593-MM05-P



Human GBA3 expression in Jurkat cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and then stained with PE Mouse anti-GBA3. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.