

CD200RLa / CD200R1L Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11620-RP02

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD200RLa / CD200R1L protein (Catalog#11620-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CD200RLa / CD200R1L (rh CD200RLa / CD200R1L; Catalog#11620-H08H; AAT00538.1; Met1-Leu239). CD200RLa / CD200R1L specific IgG was purified by Human CD200RLa / CD200R1L affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD200RLa / CD200R1L
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CD200R1L.

Please Note: Optimal concentrations/dilutions should be determined by the end user.