

Influenza A H1N1 Hemagglutinin / HA Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11684-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant H1N1 HA Protein (Catalog#11684-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/Puerto Rico/8/1934) hemagglutinin (Catalog#11684-V08H; ABD77675.1; Met 1-Gln 528). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	H1N1 (A/Puerto Rico/8/1934) HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB, ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
Western Blot	This antibody can be used at 1:500-1:1000 with the appropriate secondary reagents to detect H1N1 (A/Puerto Rico/8/1934) HA in WB.
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect H1N1 (A/Puerto Rico/8/1934) HA.

Please Note: Optimal concentrations/dilutions should be determined by the end user.