

Influenza A H10N3 Hemagglutinin / HA Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11693-RP01

GENERAL INFORMATION	
Immunogen:	11693-V08H
Preparation	Produced in rabbits immunized with purified, recombinant H10N3 (A/duck/Hong Kong/786/1979) HA / Hemagglutinin (Catalog#11693-V08H; BAF46762.1; Met1-Ile525). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	H10N3 (A/duck/Hong Kong/786/1979) HA / Hemagglutinin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	HA
APPLICATIONS	
Applications:	WB,ELISA,HemagglutininInhibition(HI)
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
Hemagglutination Inhibition	This Influenza A H10N3 (A/duck/Hong Kong/786/1979) Hemagglutinin / HA Protein neutralizing antibody can inhibit the hemagglutinin activity of recombinant Influenza A H10N3 (A/duck/Hong Kong/786/1979) Hemagglutinin / HA Protein (Cat# 11693-V08B), the HI titer is 1-20 µg/mL When 1 unit of protein added.
Western Blot	This antibody can be used at 1:500-1:1000 with the appropriate secondary reagents to detect H10N3 (A/duck/Hong Kong/786/1979) HA in WB.
ELISA	ELISA: 1:1000-1:2000

Please Note: Optimal concentrations/dilutions should be determined by the end user.