

Human Immunodeficiency Virus type 1 (HIV-1) p24 / Capsid Protein p24 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11695-R0023

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant HIV-1 Gag-P24 protein (Catalog#11695-V08E)
Preparation:	This antibody was obtained from a rabbit immunized with purified, recombinant HIV-1 Gag-P24 (rv Gag-P24; Catalog#11695-V08E; ACI05538.1; Pro 133-Leu 363).
Ig Type:	Rabbit IgG
Clone ID:	0023
Specificity:	HIV-1 (group M, subtype B, strain 92418) Gag-p24
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.