

CEACAM8 / CD66b / CD67 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11729-MM04-P

GENERAL INFORMATION

Immunogen:	Recombinant Human CEACAM8 / CD66b protein (Catalog#11729-H08H)
Reagents:	PE-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CEACAM8 / CD66b (rh CEACAM8 / CD66b; Catalog#11729-H08H; NP_001807.2; Met 1-Ser 319) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG1
Clone ID:	04
Specificity:	Human CEACAM8 / CD66b
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	CEACAM8

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CEACAM8 / CD66b / CD67 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11729-MM04-P



Profile of peripheral blood granulocytes analyzed by flow cytometry.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system.
Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>