

CALML5 Antibody (HRP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11783-MM01-H

GENERAL INFORMATION

Immunogen:	Recombinant Human CALML5 protein (Catalog#11783-H20E)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CALML5 (rh CALML5; Catalog#11783-H20E; AAH39172.1; Met 1-Glu 146). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Mouse IgG1
Clone ID:	2A3E2
Specificity:	Human CALML5
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	CALML5

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.