

CALML5 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11783-R119

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CALML5 protein (Catalog#11783-H20E)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CALML5 (rh CALML5; Catalog#11783-H20E; AAH39172.1; Met1-Glu146).
Ig Type:	Rabbit IgG
Clone ID:	119
Specificity:	Human CALML5
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:100-1:500

Please Note: Optimal concentrations/dilutions should be determined by the end user.


CALML5 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11783-R119

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of human CALML5 in human skin with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). Positive staining was localized to epidermis.