

M-CSF/CSF-1 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11792-MM16

GENERAL INFORMATION	
Immunogen:	Recombinant Human M-CSF/CSF-1 Protein (Catalog#11792-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human M-CSF/CSF-1 (rh M-CSF/CSF-1; Catalog#11792-H08H; P09603-1; Met1-Asn190). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	16
Specificity:	Human M-CSF/CSF-1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human M-CSF/CSF-1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.