

CD300C Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11832-MM03-A

GENERAL INFORMATION

Immunogen:	Recombinant Human CD300C protein (Catalog#11832-H08H)
Reagents:	APC-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD300C (rh CD300C; Catalog#11832-H08H; NP_006669.1; Met 1-Arg 183) and conjugated with APC under optimum conditions, the unreacted APC was removed.
Ig Type:	Mouse IgG1
Clone ID:	03
Specificity:	Human CD300C
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	CD300C

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD300C Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11832-MM03-A



Flow cytometric analysis of anti-human CD300c on human whole blood monocytes. The fluorescence histograms were derived from events with the forward and side light-scatter characteristics of viable monocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.