

CNTF Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11841-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human CNTF Protein (Catalog#11841-H07E)
Preparation	Produced in rabbits immunized with purified, recombinant Human CNTF (rh CNTF; Catalog#11841-H07E; NP_000605.1; Ala2-Met200). CNTF specific IgG was purified by Human CNTF affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CNTF
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	HCNTF

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CNTF.
-------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.