

Caspase-14 / CASP14 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11856-R145

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human Caspase-14 / CASP14 protein (Catalog#11856-H07E)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Caspase-14 / CASP14 (rh Caspase-14 / CASP14; Catalog#11856-H07E; NP_036246.1; Ser 2-Gln 242).
Ig Type:	Rabbit IgG
Clone ID:	145
Specificity:	Human Caspase-14 / CASP14
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,FCM
RECOMMENDED CONCENTRATION	
Flow Cytometry	FCM: 1:25-1:100
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CASP14.

Please Note: Optimal concentrations/dilutions should be determined by the end user.


Caspase-14 / CASP14 Antibody, Rabbit MAb

Catalog Number: 11856-R145



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Human CASP14 expression in HeLa cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and stained with Purified Rabbit anti-CASP14 (11856-R145), then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.